

# Understanding Your Respiratory Protection Options



Type of Respirator	CAPR	N95 Respirator	Elastomeric Half Facepiece Respirator
Testing and Approval	NIOSH approved and meets OSHA requirements		
Face Seal Fit	Loose fitting	Tight fitting	Tight fitting
Fit Testing Requirement	No	Yes	Yes
User Seal Check	No	Yes, each time the user dons (puts on) respirator	
Filtration	Filters at least 95% of airborne particles including large and small particles.		
Use limitations	Reusable (must be cleaned, disinfected and stored between each patient)	Follow current guidance.	Reusable (must be cleaned, disinfected and stored between each patient)