

PUGET SOUND REGION COVID-19 PCR—version start January 14, 2021

87798K COVID 19 PCR (MHS Labs)

Sunquest Order: COVIDI

Testing Locations: TG/MB, GSH, AMC, AH, CMC

- IP Admissions / ED Hospital Admissions
- Same Day Emergency Surgeries (<24 hours)
- Laboring OB Patients
- MHS Employee Exposures and Outbreaks: Initial testing
- Pre-Surgical <96 hours (**Pink** label w/ procedure date/time)
- Other patients by Laboratory Medical Director approval ONLY

87798S COVID 19 PCR (Batch Testing, MHS Labs)

Sunquest Order: COVIDS

Testing Location: TG/MB for all PSR

- Surveillance for MHS Employee exposure and outbreaks
- Emergency/OCED Department Dischargeable Patients
- PSR Indigo and Respiratory Clinics
- Congregate Living patients and employees
- Non-Laboring OB Patients
- Uninsured, First Responders

87798I COVID 19 PCR (UW Laboratory)

Sunquest Order: COVID

Testing Location: UW

- MHS Ambulatory patients
- Woodcreek Pediatrics
- Pre-Surgical screening (96 hrs +) (**Red** label w/ procedure date/time)

L01997 COVID 19 PCR (POCT on Roche LIAT)

Sunquest Order: N/A

Testing Location: Select POCT locations

- Available in select Indigo, WCP, OCED, and Respiratory Clinic locations for on-site testing per clinic guidelines. (Lower acuity testing collected and submitted with test code 87798S/ Batch testing.)

SPECIMEN COLLECTION Same for all orders, all locations, all laboratories

Nasal/Midturbinate (larger swab) **OR**

NP (Nasopharyngeal) (smaller swab)

(Both swabs in pack, use as indicated by provider in order and discard unused swab)

Orderable in Lawson: #1044556 (Kit CTM 2ML SARS-COV-2 NASO/OR 1KT)

**Note: Other respiratory virus tests, check for correct specimen collection instructions .

