

# PUGET SOUND REGION COVID-19 PCR Ordering and Specimen Collection Supplies

## 87798K COVID 19 PCR (MultiCare Laboratories)

Sunquest Order: COVIDM Testing Locations: TG, GSH

Use for

Most Inpatients and Higher Acuity ED patients


MHS Staff, Uninsured Patients

Ambulatory Pre-Surgical and any Pre-Surgical screening needed <48 hours before surgery

### NP Nasopharyngeal


Saline/Blue Label

FloQ swab




### Nasal/Midturbinate

Polyester swab



OR



Use one swab from pack, discard unused swab and holder

## 87798I COVID 19 PCR (UW Laboratory)

Sunquest Order: COVID Testing Location: UW

Use for


Most Ambulatory patients and Lower Acuity ED patients

First Responders

Inpatient Pre-Surgical case screening at least 48 hours prior to surgery

### Combination Swabs


2 Polyester Swabs  
Nasal and OP (throat)




Use both swabs from pack, discard unused holder

### Nasal/Midturbinate

Polyester swab



OR



Use one swab from pack, discard unused swab and holder